





```
float opt_int_time = 1.0f;          // optimal integration time, initialized to 1.0ms
int autoRepeatCounter = 0;

// This is called once to start the automatic calibration
void OnAutocalibInt()
{
    if (!g_DeviceDetected) {
        MessageBox("Device Not Attached");
        return;
    }

    int result = MessageBox("Please insert tubes for calibration.");
    if (result != IDOK) return;

    SetIntegrationTime(opt_int_time);

    StartTimer(TIMER_ID, 1000);      //Start timer to perform auto calibration
}

// Timer portion of the code
.....

        case TIMER_ID:
        {
            int done = AutocalibInt();
            if (done) KillTimer(TIMER_ID);
        }
}

.....

#define SAMPLE_READING_TARGET 600 // Target sample reading
#define REPEAT_TIMES 10

// This is called by the timer interrupt
int CGrADlg::AutocalibInt()
{
    int done = 0;

    CaptureFrame();                 // capture a frame of image

    int sample_reading;
```



The factor is inversely proportional to the initial reading:

```
if (autoRepeatCounter == 0)
{
    if (sample_reading < 10) sample_reading = 10; // avoid divide by zero
    inc_factor = 0.6f / (float)sample_reading;    // 0.6f is an empirical number
}
```

This means the smaller the initial reading, the weaker the signal, the bigger step we will take to increase integration time in proportion to sample reading difference to target level.





```
    dynIntTime = false; // clear the flag
    MessageBox ("IntTime factor adjusted to %3.2f", factorIntTime);
  }
}
```

